











Types of shanks and total lengths  
 Types de tiges et longueurs totales  
 Tipos de vástagos y longitudes totales

ID-No.	Bezeichnung	Länge mm	ø mm
 ... <b>104</b> ... ..	<b>Handstück, normal</b> Handpiece, normal length Pièce à main, longueur normale Pieza de mano, largo regular	44,00	2,35
 ... <b>105</b> ... ..	<b>Handstück, chirurgische Länge</b> Handpiece, surgical length Pièce à main, longueur chirurgicale Pieza de mano, largo quirúrgicas	65,00	2,35
 ... <b>202</b> ... ..	<b>Winkelstück, Miniatur</b> Right angle, miniature Contre angle, miniature Contrángulo, miniatura	16,00	2,35
 ... <b>204</b> ... ..	<b>Winkelstück, normal</b> Right angle, normal length Contre angle, normale Contrángulo, largo regular	22,50	2,35
 ... <b>205</b> ... ..	<b>Winkelstück, extra lang</b> Right angle, extra long Contre angle, extra-long Contrángulo, extra largo	26,50	2,35
 ... <b>206</b> ... ..	<b>Winkelstück, chirurgische Länge</b> Right angle, surgical length Contre angle, longueur chirurgicale Contrángulo, largo quirúrgicas	34,00	2,35
 ... <b>313</b> ... ..	<b>FG, shorty</b> FG, Shorty FG, Shorty FG, Shorty	16,50	1,60
 ... <b>314</b> ... ..	<b>FG, normal</b> FG, normal length FG, longueur normale FG, largo regular	19,00	1,60
 ... <b>315</b> ... ..	<b>FG, extra lang</b> FG, extra long FG, extra-long FG, extra largo	21,00	1,60
 ... <b>316</b> ... ..	<b>FG, surgical</b> FG, Surgical FG, Surgical FG, Surgical	25,40	1,60